

# Glooko data download instructions

- Can be used with many standard blood glucose monitors (such as Accu-Chek, Contour, FreeStyle, One Touch and Relion)
- You will need a cable to plug the meter into the computer. Contact your meter company at the 1-800 number on the back of your meter to request a cable.
- Create a Glooko account at [join.glooko.com](http://join.glooko.com)
- Go to [my.glooko.com](http://my.glooko.com), download the Glooko Uploader installation file, download the meter
- To share data with your provider:
  - Enter our ProConnect Code (childrensdayton) during set-up
  - Notify office staff anytime you upload, since a notification **is not** sent to your provider.

For more information and the most current instructions, ask us or visit [my.glooko.com](http://my.glooko.com).

# Dexcom data download instructions

- If using a phone as a receiver, install the Dexcom Clarity app on your phone.
- If using a receiver, download the receiver to the Dexcom Clarity webpage ([clarity.dexcom.com](http://clarity.dexcom.com)).
- To share data with your provider, you can:
  - Get office share code from office staff **or**
  - Generate 12 letter share code, then give that code to office staff.
  - Notify office staff anytime you upload, since a notification **is not** sent to your provider.

For more information and the most current instructions, ask us or visit [clarity.dexcom.com](http://clarity.dexcom.com)

# FreeStyle Libre data download instructions

- Create a LibreView account at [libreview.com](http://libreview.com) or through the Freestyle LibreLink app.
  - If using a phone as the reader, install the LibreLink app on your phone.
  - If using a reader, download the reader to the LibreView webpage ([libreview.com](http://libreview.com)).
- To share data with your provider, you can:
  - Enter Practice ID (13098337) during set-up.

OR

- Give your email to office staff, accept email invitation to share your data with your provider.
- Notify office staff any time you upload since a notification **is not** sent to your provider.

For more information and the most current instructions, ask us or visit [pat.libreview.io](http://pat.libreview.io).

# Medtronic data download instructions

- Create a CareLink account at [carelink.minimed.com](http://carelink.minimed.com).
- Download Carelink Uploader installer.
- Log in to Carelink, click “Upload Device” and follow the instructions.
- To share data with your provider:
  - Give your provider your CareLink username and password to allow them to access your CareLink report.
  - Notify office staff any time you upload since a notification **is not** sent to your provider.

For more information and the most current instructions, ask us or visit [medtronicdiabetes.com/products/carelink-personal-diabetes-software](http://medtronicdiabetes.com/products/carelink-personal-diabetes-software).

# Omnipod data download instructions

- Eros pump: Go to [join.glooko.com](http://join.glooko.com) to create an account, then [my.glooko.com](http://my.glooko.com) to upload the pump. The pump PDM must be plugged into the computer to download.
- DASH pump: Download the Glooko app, then open the app and click “Sign up for Glooko”. At setup, do the following so the PDM can auto upload every 24 hours:
  - Create Podder Central Account with Omnipod.
  - Create Glooko account.
  - Make sure wifi is on.
  - Enable auto uploads with Glooko web browser/Link Podder Central.
  - Enable auto upload via Glooko app.
- To share data with your provider:
  - Enter our ProConnect Code (childrensdayton) during set-up.
  - Notify office staff any time you upload since a notification **is not** sent to your provider.

For more information and the most current instructions, ask us or visit [join.glooko.com](http://join.glooko.com).

# Tandem data download instructions

- Create t:connect account at [tconnect.tandemdiabetes.com](http://tconnect.tandemdiabetes.com).
- On your first time uploading, plug the pump into the computer and upload on the t:connect website
- After that, the pump can upload automatically using the t:connect mobile app. Turn on your Bluetooth settings, and connect to wifi.
- To share data with your provider:
  - Give provider your t:connect account email and password to allow them access to your t:connect report.
  - Notify office staff any time you upload since a notification **is not** sent to your provider.

For more information and the most current instructions, ask us or visit [tandemdiabetes.com/products/software-apps/tconnect-application](http://tandemdiabetes.com/products/software-apps/tconnect-application).

# One Touch Verio Flex/Reflex data download instructions

- Download the One Touch Reveal Mobile App
- Download the meter results remotely to the app.
- To share data with your provider:
  - Send the logbook PDF via MyKidsChart **or**
  - Enter the clinic code, HHQFXW
  - Notify office staff any time you upload since a notification **is not** sent to your provider.

For more information and the most current instructions, ask us or visit <https://www.onetouch.com/products/glucose-meters/onetouch-verio-flex>