

# Omnipod data download instructions

- Eros pump: Go to [join.glooko.com](http://join.glooko.com) to create an account, then [my.glooko.com](http://my.glooko.com) to upload the pump. The pump PDM must be plugged into the computer to download.
- DASH pump: Download the Glooko app, then open the app and click “Sign up for Glooko”. At setup, do the following so the PDM can auto upload every 24 hours:
  - Create Podder Central Account with Omnipod.
  - Create Glooko account.
  - Make sure wifi is on.
  - Enable auto uploads with Glooko web browser/Link Podder Central.
  - Enable auto upload via Glooko app.
- To share data with your provider:
  - Enter our ProConnect Code (childrensdayton) during set-up.
  - Notify office staff any time you upload since a notification **is not** sent to your provider.

For more information and the most current instructions, ask us or visit [join.glooko.com](http://join.glooko.com).