

donning and doffing

droplet and contact isolation for COVID-19





droplet/contact + face shield

- This isolation is for COVID rule out or COVID positive patients without aerosolized therapies.
- This isolation type will require the following PPE (personal protective equipment):
 - Gown
 - Surgical mask
 - Gloves
 - Face shield

donning

1. Perform hand hygiene	
2. Put on isolation gown <ul style="list-style-type: none">• Ensure that gown is secured at back of neck• Tie gown at waist and ensure no gaps in the back	
3. Apply surgical mask <ul style="list-style-type: none">• The same surgical mask will be worn for your entire shift• If soiled, change immediately	
4. Put on face shield	
5. Apply gloves	
6. Ready for patient care	

doffing droplet/contact + field shield

<p>1. Remove dirty gown INSIDE patient room</p> <ul style="list-style-type: none">• Take care to only touch outside of the dirty gown with your gloves.	
<p>2. Remove dirty gloves INSIDE patient room</p>	
<p>3. Sanitize your hands</p>	
<p>4. Proceed to cleaning station with face shield and mask still on</p>	
<p>5. Cleaning face shield</p> <ul style="list-style-type: none">• Perform hand hygiene• Don new gloves• Obtain one Oxivir wipe from outside cabinet• Remove face shield (see picture to right)• Clean inside of face shield• Clean outside of face shield• Wipe off table using same Oxivir wipe• Discard wipe• If removing mask, remove mask at this point in process	
<p style="text-align: center;">DO NOT DISCARD MASK IF NOT SOILED *For Oxivir to be effective, item must remain wet for 60 seconds. After 60 seconds, item may be wiped dry if needed. There is no mandated dry time.</p>	